

## E.Z.N.A. Fungal DNA Extraction Mini-Prep.

1. **Add a pinch of sand** to each tube along with **200 µl of buffer FG1**. Let sit for 5 min then another 10 min at **-20°C to soften the tissue**.
2. **Remove tubes from -20°C**, let thaw for 5 min, then **grind the tissue with a micropestle** until all chunks/pieces are homogenized.
3. **Add 600 µl buffer FG1**.
4. **65°C for 10 min**. Mix 2X.
5. **180 µl buffer FG2**. Vortex.
6. 5 min on ice.
7. **Centrifuge 10 min @ 10K g**. BALANCE and HINGES UP ALWAYS!
8. **Transfer 700 µl / 600 µl supernatant to new tube**. Do not transfer any tissue!
9. **0.7 vol (480 µl / 420 µl) isopropanol**. Vortex.
10. **Centrifuge 2 min @ 10K g**.
11. **Remove supernatant**. Do not disturb DNA pellet (if visible). Invert on paper towel to let dry for 1 min.
12. **300 µl ddH<sub>2</sub>O heated to 65°C**. Vortex.
13. **4 µl RNase A**. Vortex.
14. **150 µl buffer FG3 + 300 µl 100% ETOH**. Vortex.
15. **Transfer** to HiBind Mini Column and supplied collection tube.
16. **Centrifuge 1 min @ 10K g**. Discard filtrate.
17. 750 µl DNA Wash Buffer.
18. **Centrifuge 1 min @ 10K g**. Discard filtrate.
19. Repeat steps 17 and 18.
20. **Centrifuge 2 min @ 10K g**. To remove all liquid.
21. **Transfer Mini Column to new 1.5µl microcentrifuge collection tube**. Let sit for 5-10 minutes so the wash buffer can evaporate.
22. **50 µl Elution Buffer heated to 65°C**. (Optional: let sit for 3-5 minutes).
23. **Centrifuge 1 min @ 10K g**.
24. Repeat step 22-23 one more time.

## E.Z.N.A. Fungal DNA Extraction Mini-Prep.

1. **Add a pinch of sand** to each tube along with **200 µl of buffer FG1**. Let sit for 5 minutes then another 10 at **-20°C to soften the tissue**.
2. **Remove tubes from -20°C**, let thaw for 5 min, then **grind the tissue with micropestle** until all chunks/pieces are homogenized.
3. **Add 600 µl buffer FG1**.
4. **65°C for 10 min**. Mix 2X.
5. **180 µl buffer FG2**. Vortex.
6. 5 min on ice.
7. **Centrifuge 10 min @ 10K g**. BALANCE and HINGES UP ALWAYS!
8. **Transfer 700 µl / 600 µl supernatant to new tube**. Do not transfer any tissue!
9. **0.7 vol (480 µl / 420 µl) isopropanol**. Vortex.
10. **Centrifuge 2 min @ 10K g**.
11. **Remove supernatant**. Do not disturb pellet. Invert on paper towel to let dry for 1 min.
12. **300 µl ddH<sub>2</sub>O heated to 65°C**. Vortex.
13. **4 µl RNase A**. Vortex.
14. **150 µl buffer FG3 + 300 µl 100% ETOH**. Vortex.
15. **Transfer** to HiBind Mini Column and supplied collection tube.
16. **Centrifuge 1 min @ 10K g**. Discard filtrate.
17. 750 µl DNA Wash Buffer.
18. **Centrifuge 1 min @ 10K g**. Discard filtrate.
19. Repeat steps 17 and 18.
20. **Centrifuge 2 min @ 10K g**. To remove all liquid.
21. **Transfer Mini Column to new 1.5µl microcentrifuge collection tube**. Let sit for 5-10 minutes so the wash buffer can evaporate.
22. **50 µl Elution Buffer heated to 65°C**. (Optional: let sit for 3-5 minutes).
23. **Centrifuge 1 min @ 10K g**.
24. Repeat step 22-23 one more time.