

## QuickNDirty Extraction Solutions

Note: This protocol has been handed down from lab to lab. The Wilson Lab received this version, with protocols adjusted to fit our lab, from Dr. Brian Perry. October 2019.

### Extraction Solution (ES)

Stocks: 1 M Tris solution (pH=8.0)

- 1) add 10 ml of 1 M Tris stock into clean 100 ml vessel
- 2) add 1.86 g KCl
- 3) add 0.37 g EDTA
- 4) add 80 ml DI or ultrapure H<sub>2</sub>O and shake until solutes dissolve
- 5) Titrate with 1 M NaOH to pH = ~ 9.5-10.0
- 6) top up to 100 ml with DI H<sub>2</sub>O
- 7) filter sterilize and aliquot into sterile 1.7 mL eppendorf tubes (filtering is recommended but optional)

### Dilution Solution (BSA 3%)

- 1) add 3 g of BSA (e.g. Sigma or Omni – 98-99% purity, heat shock fractionated) into clean vessel
- 2) top up to 100 ml with DI/Ultrapure H<sub>2</sub>O
- 3) shake BSA into solution
- 4) filter sterilize into sterilized 2 ml eppies (same as above here)

## Procedure

1. Take a **small piece** of fungal material, not greater than one kernel of long-grain rice, and place that into an appropriately labeled 1.5 mL microcentrifuge tube.
2. Add a small pinch of **sterile sand** ([Fisher Sci S25-500](#)) to each tube.
3. Pipette out **100µL of ES** into each tube. Let tissue absorb ES for 5 minutes.
4. Place tubes into -20°C **freezer** until tissue is frozen (~15 min). This softens the tissue for grinding. Remove and thaw.
5. **Grind** with sterilized micropestle ([Fisher Sci 12-141-368](#)) until all visible chunks of tissue are broken down and homogenized.
6. **Incubate** at room temp for 10+ minutes then incubate for 10 minutes at 95° C.
7. Add an equal volume of **Dilution Solution (DS)** (3% BSA).
8. **Vortex** and **Spin** briefly (max for 30 sec).
9. Samples are now ready for PCR. 1-2 µL is usually sufficient for PCR (from a 20µL ES + 20µL DS = 40µL volume extraction). For herbarium specimens we always do a 1/10 to 1/1000 dilution series, and use 2 µL of each for PCR.
10. Store DNA extractions in freezer.

NOTE: It's currently not clear how long these DNA extractions will last in the freezer for. Current discussions have suggested that after 6 months, DNA that was extracted using this method breaks down is no longer viable. Further exploration and testing is needed.